Reorder Number

REF

5551000

REF

9551000

Refurbished

Non-sterile

X X-STREAM™

Laparoscopic Irrigation Controller

Operator's Manual**Guide de l'utilisateur****Caution:** Read all instructions prior to use.**BAIRD**

DAVOL INC.

PK3797715
142R

	Safety Information and Symbols Definition	en 1
	Accessories	en 1
Section 1	Introduction	en 2
Section 2	Warnings	en 2
	A. X-STREAM™ Laparoscopic Irrigation Controller	
	B. X-STREAM™ Laparoscopic Irrigation Tubing Set	
Section 3	Precautions	en 2
Section 4	Specifications	en 3
Section 5	User Interface (Description/Function)	en 3
Section 6	Controller Setup and Operating Instructions	en 4-5
Section 7	Calibration/Verification Testing	en 5
	A. Grounding Resistance	
	B. System Calibration Procedure	
	C. Conclusions	
Section 8	Troubleshooting	en 6
Section 9	Cleaning and Maintenance	en 6
Section 10	Controlling System Audio Indicator	en 7
Section 11	Appendix	en 8
	A. Return Policy Notification	
	1. Authorization	
	2. Freight (Returned Goods)	
	3. Merchandise Returnable for Credit	
	B. Limited Product Warranty	
	C. Product Changes	
	D. Excuses	
	E. Safety Agency Testing	
	Product Illustration	en 9

SAFETY INFORMATION and SYMBOLS DEFINITION

CAUTION: Read all instructions prior to use.

The exclamation point within an equilateral triangle is intended to alert the user to the presence of important operating and maintenance (servicing) instructions in the literature accompanying the device.

The lightning flash with arrowhead symbol, within an equilateral triangle, is intended to alert the user of the presence of uninsulated "dangerous voltage" within the product's enclosure that may be of sufficient magnitude to constitute a risk of electric shock to persons.

Maximum flow rate setting

Other Patents Pending

Type B Applied Part

Default flow rate setting

Fuse

IPX4

Protected against the effects of splashing water.

PRIME Prime

Reusable

Warranty void if seal is broken.

"WARNING": Risk of fire.
Replace fuse as marked.

Contents

Refurbished

Indicates separate collection for electrical and electronic equipment.

Biological Risks

Temperature Limitation

Consult Instructions for Use

Protective earth (ground)

Non-Sterile

Contains phthalates

ACCESSORIES

REF 5552000

X-STREAM™ Laparoscopic Irrigation Tubing Set with SMOKEVAC™ Trumpet Valve, 5mm x 33cm Probe Tip with Holes and Cojoined Suction/Irrigation Tubing

REF 5552002

X-STREAM™ Laparoscopic Irrigation Tubing Set with NEZHAT-DORSEY™ Trumpet Valve, 5mm x 33cm Probe Tip with Holes and Cojoined Suction/Irrigation Tubing

REF 5552001

X-STREAM™ Laparoscopic Irrigation Tubing Set with SMOKEVAC™ Trumpet Valve and Cojoined Suction/Irrigation Tubing, without Probe Tip

REF 5552003

X-STREAM™ Laparoscopic Irrigation Tubing Set with NEZHAT-DORSEY™ Trumpet Valve and Cojoined Suction/Irrigation Tubing, without Probe Tip

The X-STREAM™ Laparoscopic Irrigation Controller is a versatile and reusable electromechanical controller (Figure 1, Page en 3). It is small, lightweight, and pole-mounted. It has two easily adjustable flow rate settings, and may be used with any of the DAVOL™ X-STREAM™ Laparoscopic Irrigation Tubing Sets.

SECTION 2 WARNINGS

A. X-STREAM™ Laparoscopic Irrigation Controller

1. The performance of the X-STREAM™ Laparoscopic Irrigation Controller depends on proper setup. Please refer to tubing set Instructions for Use for setup instructions.
2. Use of any modified or unauthorized disposables with the X-STREAM™ Laparoscopic Irrigation Controller is contraindicated and may result in undesirable system performance.
3. If performing electrical safety testing follow IEC 60601-1:2006 recommendations.
4.  The X-STREAM™ Laparoscopic Irrigation Controller contains internal components that could cause electrical shock hazard. Do not open the housing of the Controller or attempt to service the system.
5. The power supply cord is not rewirable. Please return X-STREAM™ Laparoscopic Irrigation Controller to DAVOL™ Inc. for servicing.
6. To avoid the possibility of electrical shock, always remove the tubing set and disconnect the power cord prior to cleaning.
7. Do not use the X-STREAM™ Laparoscopic Irrigation Controller in the presence of a FLAMMABLE ANESTHETIC MIXTURE WITH AIR or with OXYGEN OR NITROUS OXIDE.
8. Do not autoclave or steam sterilize the X-STREAM™ Laparoscopic Irrigation Controller or subject the pump to temperatures in excess of 70°C. To avoid possibility of electric shock, always disconnect the power cord prior to cleaning.
9. Do not touch the X-STREAM™ Laparoscopic Irrigation Controller and the patient simultaneously.
10. No user serviceable items. Return X-STREAM™ Laparoscopic Irrigation Controller to DAVOL™ Inc. for all repairs.
11. The tubing set will provide gravity flow if the pumping chamber is not properly loaded into the controller's pumping chamber receptacle.
12. Do not modify this equipment without authorization of DAVOL™, as it may compromise the safety and effectiveness of the device.
13. If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe and effective use of the equipment.
14. To ensure safety and accuracy, contact DAVOL™ for instructions regarding the correct replacement of interchangeable or detachable parts.
15. X-STREAM™ should be installed in accordance with the EMC information provided in this document.
16. If X-STREAM™ is installed in the vicinity (see tables) of a RF communication device proper operation of the X-STREAM™ should be observed.
17. Use of accessories other than declared by DAVOL™ may negatively impact X-STREAM's™ emissions and immunity.

B. X-STREAM™ Laparoscopic Irrigation Tubing Set

1. As with any surgical procedure where extensive irrigation is utilized, intravasation levels should be closely monitored by the attending physician.
2. The pressure accuracy of the X-STREAM™ System depends on proper setup. It is important that the bottom of the X-STREAM™ Controller be positioned at least 12 inches above the highest height of the patient to provide sufficient pressure for optimal performance of the system.
3. The X-STREAM™ Irrigation System is capable of generating high irrigation impact forces and flow rates. Use the appropriate, required pressure necessary at all times to accomplish the procedure.
4. The X-STREAM™ Irrigation System will generate even higher irrigation impact forces (compared to impact forces generated via a 5mm standard probe tip) when used with the DAVOL™ 10mm ENDO-POOL™, 10mm Hydro-Dissection Channel Probe or the 5mm Micro Probe tips.
5. Electrosurgery instruments and electrode inserts must not be used with uninsulated metallic probe tips. Electrosurgery instruments and inserts should be used only with available DAVOL™ non-conductive or insulated probe tips.
6. When using electrosurgical attachments or inserts, follow all instructions, warnings, precautions, and contraindications provided with these devices.

SECTION 3 PRECAUTIONS

1. Read all the instructions provided with the X-STREAM™ Laparoscopic Irrigation Tubing Set prior to use.
2. The controller automatically turns "ON" when the tubing set is loaded into the controller. The logo and "PRIME" LED's will light to indicate that the controller is on. There is no ON/OFF switch.
3. Do not remove the square, metal contact label located at the top back of the tubing set pumping chamber.
4. Solvents should not be used on the controller's logo and keypad as they could cause damage to the device.
5. When air enters the system, the controller returns to the "PRIME" state, the "PRIME" text will blink, and the audio indicator will sound.
6. Gravity prime and perform all priming with the probe tip outside of patient.
7. Please use special precautions regarding the use of the X-STREAM™, The X-STREAM™ is to be used in EMC environments that are suitable for Group 1, Class A medical equipment.
8. Grounding reliability can only be achieved when power supply is connected to an equivalent receptacle marked "Hospital Only" or "Hospital Grade". [US only]
9. X-STREAM™ should be positioned so that it is not difficult to disconnect device from the power source, in order to facilitate access.
10. The X-STREAM™ controller should not be placed adjacent to or stacked with other equipment unless nominal operation of the X-STREAM™ has been observed as it may cause a malfunction of the device.

Physical	Electrical	System	Environmental/Storage
Width: 16.0cm	Voltage: 100-240VAC	Classification: Class 1	Generated Noise: 85dBA
Length: 15.2cm	Frequency: 50/60 Hz	Mode: Intermittent Operation	Ambient Temperature: 15°C to 30°C
Height: 18.4cm	Fuse(s): 2.5A external	Type:  Type B Applied Part	Relative Humidity: 30% - 65% (Non-Condensing)
Weight: 1.2kg		Ingress Protection: IPX4	Atmospheric Pressure: 70kPa to 100kPa
		Continuous Irrigation Time (max): 15 sec/cycle	

SECTION 5 USER INTERFACE

ITEM	DESCRIPTION	FUNCTION
1	DAVOL logo LED	An illuminated logo indicates the presence of a properly loaded disposable tubing set and that the controller is "ON", waiting for "PRIME" or "READY".
2	PRIME LED	Illuminated "PRIME" text indicates that the tubing set needs to be primed.
3	"HI" setting LED button	Depressing this button activates the maximum flow rate. There are only two flow rate settings.
4	"LO" setting LED button	Depressing this button activates the default flow rate. There are only two flow rate settings.
5	IV pole clamp	Mounts controller to standard IV poles.
6	Fuse holder	Holds two 2.5 Amp. fuses. Lift tab on inside of IEC outlet to access fuse holder.
7	IEC power cord	Attachable power cord compatible with IEC power cord receptacle. Attaches to (8).
8	IEC power cord receptacle	Female receptacle for IEC compatible, attachable power cords.
9	Power cord retention strap stud	For securing the strap.
10	Power cord retention strap	Useful in retaining the power cord during storage of the controller.
11	Metal case	 Grounded surface for testing leakage current.
12	Pumping chamber receptacle	Receptacle to receive disposable pumping chamber.

Figure 1



Figure 2

Important: Thoroughly review all of the operating instructions in this manual and in the X-STREAM™ Laparoscopic Irrigation Tubing Set Instructions for Use before using the X-STREAM™ Laparoscopic Irrigation Controller.

X-STREAM™ Controller Setup

1. With the bottom of the X-STREAM™ Controller at least 12 inches above the highest height of the patient at any time during surgery, mount the controller on the IV pole. To mount, loosen pole clamp knob by turning it. Position the pole clamp around the IV pole and tighten the knob to secure the controller to the pole.
2. Ensure a IEC/Hospital Grade power cord is plugged into the IEC receptacle on the back of the controller. Also ensure that the power cord is connected to an operable Hospital Grade AC outlet.
3. Hang irrigation bag(s) on the IV pole. Because of the high flow rates generated by the X-STREAM™ Irrigation System, 2000cc or larger IV bag sizes are recommended to minimize bag changes.
4. If spiking two bags is desirable, use the optional Dual Spike Adapter (REF 0026110).

Note: Changing the setting is possible both with and without the motor running.

Note: The flow rate setting cannot be changed when the controller is in the “PRIME” mode.

Precaution: Do not remove the square, metal contact label located at the top back of the tubing set pumping chamber.

Opening the Tubing Set Package

1. The circulating nurse opens the package and delivers the tubing set onto the sterile field.
2. The scrub nurse retains the trumpet valve handpiece and checks and tightens the knurled Quick Disconnect Adapter at the front of the Trumpet Valve. The scrub nurse securely attaches a probe tip by pushing the probe tip over the Quick Disconnect Adapter (on some tips a double “click” is heard), then hands the pumping chamber (X-STREAM™ clear cartridge) to the circulating nurse.
3. The circulating nurse connects the suction tubing to the suction source. Peel away from irrigation tubing as needed.

Loading the Tubing Set

Precaution: Read all directions for the X-STREAM™ Laparoscopic Irrigation Tubing Set before beginning the loading procedure.

1. Slide the pumping chamber downward into the controller until it snaps into place flush with the top and face of the controller.
Precaution: The controller automatically turns “ON” when the tubing set is loaded into the controller. The logo and “PRIME” LED’s will light to indicate that the controller is on. There is no ON/OFF switch.
2. Remove spike cap and spike the irrigation bag. Leave sufficient slack in line to facilitate bag placement/removal. If spiking two bags is desirable, use the optional Dual Spike Adapter (REF 0026110).
3. Refer to X-STREAM™ Laparoscopic Irrigation Tubing Set Instructions for Use for information on the probe handle and tips.

Priming the System

Note: The system includes audible and visual indicators to let the user know priming is required. The “PRIME” LED light will flash and the audible indicator will sound when the system detects the presence of the pumping chamber upon initial setup. See Section 10, page en 7, for controlling the audio indicator.

Precaution: Gravity prime and perform all priming with probe tip outside of patient.

1. To prime, depress irrigation button (colored button) until all air is out of the system and only irrigant exits the probe tip.
2. Release the irrigation button.
3. If the pumping chamber is not primed properly, the “PRIME” text will continue to flash and the audio indicator will continue to sound after the irrigation button is released.

Note: Removal and reloading of the pumping chamber will return the controller to the “PRIME” mode. The DAVOL™ logo LED will illuminate to indicate that the X-STREAM™ Controller is “ON”. The “PRIME” text will flash and the audio indicator will sound to indicate that the tubing set needs to be primed.

4. When ready to irrigate, depress the irrigation button (colored button) on the trumpet valve.

Note: The motor will run when the irrigation button is depressed and will stop after it is released. A safety feature turns the pump off after 15 seconds of continuous use. Release the irrigation button and then re-depress to reactivate the pump.

Basic Operation

1. Changing the Flow Rate Setting

1. There are two flow rate settings available on the X-STREAM™ Laparoscopic Irrigation Controller. The LED lights to identify the activated flow rate setting.
2. The controller automatically activates at the low setting. To select the maximum flow rate setting, the circulating nurse should press the “HI” button.
3. Switch back and forth between the “LO” and “HI” settings throughout the procedure by pressing the appropriate button.
4. The surgeon can control the flow at the handpiece.

**Default
Flow Rate**
2500 - 3800 cc/min

**Maximum
Flow Rate**
3400 - 5000 cc/min

Note: Changing the setting is possible both with and without the motor running.

Note: The flow rate setting cannot be changed when the controller is in the “PRIME” mode.

Note: The NEZHAT-DORSEY™ Trumpet Valves offer variable flow control. For example, depressing the irrigation button on the handpiece halfway delivers approximately one-half of the available irrigation pressure/flow.

2. Changing Irrigation Bags

Note: The system includes audible and visual indicators to let the user know priming is required. The “PRIME” and “HI” LED lights will blink and the audio indicator will sound when the bag runs out of fluid.

Precaution: Perform all priming with the probe tip outside the patient. If the bag(s) runs out of irrigant, the motor will stop running. When air enters the system, the controller returns to the “PRIME” state, the “PRIME” and “HI” text will blink, and the audio indicator will sound.

- 1. Hang a new bag(s) and spike.
- 2. Depress the irrigation button (colored button) on the handpiece until fluid fills the pumping chamber, and then, while continuing to hold the irrigation button down, press the “HI” button on the controller to power prime the system. Release irrigation button after 3 seconds.

Note: The “PRIME” text and the audio indicator will go off and the flow rate setting LED will light up. The system returns to the flow rate setting in use prior to bag runoff.

If the “PRIME” text remains lit, and the audio indicator is sounding, the tubing set has not been primed properly. Continue power priming to remove all air from the system.

3. Removing the Tubing Set

- 1. If using a dual spike adapter (**REF 0026110**) close the appropriate clamp(s) on the bag spike arms of the disposable tubing set.
- 2. Remove suction tube from suction device.
- 3. Remove the pumping chamber from the controller by sliding it upwards.
- 4. Discard the disposable tubing set in accordance with hospital waste disposal practices.

 **Precaution:** After use, the tubing set, trumpet valve and probe tip are a potential biohazard. Handle and dispose of as required by hospital policy and applicable laws.

- 5. Wipe surface of controller with mild soap and water solution applied with a damp cloth. Store for future use.

SECTION 7 CALIBRATION / VERIFICATION TESTING

A. Grounding Resistance

Power cord integrity may be checked without opening the controller casing by measuring grounding resistance from the Power Entry Module Metal Casing on the back of the controller to the ground pin of the attachment plug per NFPA 99 Standard for Health Care Facilities (1996 edition) Chapter 7-5.1.3.2 Resistance. All other surfaces are isolated from ground. Perform electrical safety testing following IEC 60601-1:2006 recommendations.

B. System Calibration Procedure

Electrical Voltage Verification Protocol

Purpose


To verify that the reusable controller is delivering appropriate voltage to the disposable tubing set.

Equipment

RMS Multimeter

Procedure

- 1. Set up disposable tubing set and prime system according to Instructions for Use.
- 2. Maintain a full 3L IV bag at the height described in the Instructions for Use for the disposable tubing set.
- 3. On the RMS Multimeter, connect conductor leads to Common (-) and Volts/Ohms (+).
- 4. Set dial on RMS Multimeter to measure DC Volts (\overline{V})
- 5. Contact the probed end of the conductor leads to the metal contacts at the bottom of the pumping chamber section of the disposable tubing set as shown in the **Electrical Test Setup Diagram**.
- 6. Activate the **LO** button on the controller. Irrigate and record maximum voltage reading on RMS Multimeter while irrigating.
- 7. Activate the **HI** button on the controller. Irrigate and record maximum voltage reading on RMS Multimeter while irrigating.

Calibration Table

Controller Setting	Acceptable Range of DC Voltage (± 1VDC)
DEFAULT	19 VDC
MAXIMUM	23 VDC

Electrical Test Setup Diagram

C. Conclusions

Any deviations from the test results presented above should be reported to the customer service department of the manufacturer, Davol Inc., 100 Crossings Boulevard, Warwick, RI 02886, 1-800-556-6275. If outside of the U.S.A., deviations should be reported to your local X-STREAM™ System distributor.

PROBLEM	POSSIBLE SOLUTIONS / ACTIONS
Logo LED's do not illuminate when disposable pumping chamber is installed into pumping chamber receptacle.	<ul style="list-style-type: none"> • Connect power cord securely to the IEC power cord receptacle on the back of the controller. • Connect power cord securely to an operable Hospital Grade AC outlet. • Load the pumping chamber of the disposable tubing set properly into the pumping chamber receptacle of the controller (Press firmly downward until the disposable is flush with the top surface of the controller.) • Check fuses and contact DAVOL™ for instructions regarding the correct replacement of interchangeable or detachable parts. • Check the back of the pump chamber to ensure that metallic label is installed. • Install new disposable unit. • Call DAVOL™ Customer Service if problem persists.
Logo LED's blink after disposable pumping chamber is installed and pump will not start.	<ul style="list-style-type: none"> • System fault detected. Electrical continuity to the disposable pumping chamber may be interrupted. • Install new disposable unit. • Call DAVOL™ Customer Service if problem persists.
"PRIME" text blinking after system is gravity primed.	<ul style="list-style-type: none"> • Depress irrigation button on trumpet valve. Prime unit to remove all trapped air. • Call DAVOL™ Customer Service if problem persists.
Logo LED's are illuminated but fluid is not flowing through the system.	<ul style="list-style-type: none"> • Check that the irrigation bag(s) are not empty. • Check that the irrigation bag(s) are spiked completely. • Open the clamp(s) for the appropriate irrigation bags. • Remove any kinks in the flexible tubing. Ensure that the fluid path is open and fluid is flowing toward the trumpet valve handpiece. • Check that the appropriate button (colored button) is depressed on the trumpet valve handpiece.
Low pressure and flows.	<ul style="list-style-type: none"> • Ensure that the pumping chamber is properly seated into the pumping chamber receptacle of the controller. • Ensure that spike is inserted fully into the irrigation container. • Completely open the clamps on the irrigation bag spike arms. • Remove any kinks in the flexible tubing. Ensure that the fluid path is open and fluid is flowing toward the trumpet valve handpiece. • Install new disposable unit. • Call DAVOL™ Customer Service if problem persists.
Irrigation pressure does not increase or decrease when pressure settings are selected from the keypad.	<ul style="list-style-type: none"> • Ensure that the pumping chamber is properly seated into the pumping chamber receptacle of the controller. • Completely open the clamps on the irrigation bag spike arms. • Remove any kinks in the flexible tubing. Ensure that the fluid path is open and fluid is flowing toward the trumpet valve handpiece. • Firmly depress irrigation button. • Call DAVOL™ Customer Service if problem persists.
Motor does not run. Logo and pressure setting buttons are illuminated.	<ul style="list-style-type: none"> • Check that the irrigation bag(s) are not empty. • Remove any kinks in the flexible tubing. Ensure that the fluid path is open and fluid is flowing toward the trumpet valve handpiece. • Ensure that the irrigation button on trumpet valve handpiece is depressed. • Install new disposable unit.
Motor runs intermittently without cause.	<ul style="list-style-type: none"> • Depress irrigation button. Prime unit to remove all trapped air. • Avoid violent shaking of the irrigation tubing during use.
Motor runs continuously without cause.	<ul style="list-style-type: none"> • Depress irrigation button. Prime unit to remove all trapped air. • Check for fluid leaking at the trumpet valve handpiece or irrigation tip. • Check for fluid leaking at the pumping chamber. • Call DAVOL™ Customer Service if problem persists.

SECTION 9 CLEANING AND MAINTENANCE

1. The X-STREAM™ Laparoscopic Irrigation Controller may be cleaned with a mild soap and water solution applied with a damp cloth.
2. If contacts become coated with deposits from irrigation fluid, they must be cleaned using a damp swab.
WARNING: Do not autoclave or steam sterilize the X-STREAM™ Laparoscopic Irrigation Controller or subject the controller to temperatures in excess of 70°C. To avoid possibility of electric shock, always disconnect the power cord prior to cleaning.
3. The line cord is the supply (MAINS) disconnect (ISOLATION) device.

A. Return Policy Notification

If the controller does not function after troubleshooting tips, or to have DAVOL™ service the controller, please call DAVOL™ Inc. at 1-800-556-6275 for return goods authorization.

Do not attempt to service the controller. Servicing the controller voids the warranty of the device. The controller contains a seal that will be broken if the device is opened. If the seal is broken when the controller is returned to DAVOL™, the warranty is void.

1. Authorization

Unauthorized returns will not be accepted and will be returned to the customer at the customer's expense.

A. Request for Return Authorization

The request for Return Authorization form must include:

1. the catalog number and quantity,
2. the reason for return (defective merchandise, DAVOL™ shipping error, or other),
3. the customer purchase order number and date, and
4. the DAVOL™ invoice number and date.

Note: The Return Goods Authorization will expire 180 days from the date of issue.

2. Freight (Returned Goods)

When authorized, all merchandise must be returned FREIGHT PREPAID. Any merchandise returned freight collect will be refused by DAVOL™ and returned to the customer at the customer's expense. DAVOL™ will credit freight charges to the customer for the following returns:

1. defective product, and
2. merchandise shipped in error by DAVOL™.

3. Merchandise Returnable for Credit

Merchandise will be considered for return only in minimum shipping quantities. Only items in the current price list will be considered for credit. All credit will be issued based on original net purchase price; or, if invoice number is not provided, credit will be issued based on the previous list at best price. Merchandise can be returned up to 15% of annual sales.

B. Limited Product Warranty

DAVOL™ Inc. ("Daval") warrants to the original purchaser that this product will be free from defects in material and workmanship for a period of one (1) year from the date of purchase. If this product proves to be so defective, purchaser may return same to DAVOL™ for repair or replacement, at DAVOL'S™ option. All returns must be authorized in advance in accordance with DAVOL'S™ Returned Goods Policy found in its then current Price List. The liability of DAVOL™ under this limited product warranty does not extend to any abuse or misuse of this product or its repair by anyone other than an authorized DAVOL™ representative.

[THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED (INCLUDING, WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE). THE LIABILITY AND REMEDY STATED IN THIS LIMITED PRODUCT WARRANTY WILL BE THE SOLE LIABILITY OF DAVOL™ AND REMEDY AVAILABLE TO PURCHASER FOR THIS PRODUCT, WHETHER IN CONTRACT, TORT (INCLUDING NEGLIGENCE) OR OTHERWISE, AND DAVOL WILL NOT BE LIABLE TO PURCHASERS FOR ANY INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF ITS HANDLING OR USE].

C. Product Changes

All products and product specifications identified in the Price List are based upon information available to DAVOL™ Inc., at the time of publication. The right is reserved to discontinue any of these products or to change any such specifications with or without prior notice.

D. Excuses

DAVOL™ Inc., will not be liable for its failure to perform or for a delay in performance of any order due to strikes, fire, explosion, flood, riot, lock out, injunction, interruption of transportation, unavoidable accidents, acts of government or a public enemy, inability to obtain supplies at reasonable prices or other causes beyond its control.

None of the foregoing terms may be modified except upon the express written consent of DAVOL™ Inc., stating that it is an amendment or modification thereto.


E. Safety Agency Testing

The X-STREAM™ Laparoscopic Irrigation Controller and electrical accessories have been inspected, tested and certified and are in compliance with the applicable requirements of the Standards for Medical Electrical Equipment, IEC 60601-1, UL 60601-1 and CSA C22.2 NO. 601.1.

The X-STREAM™ Laparoscopic Irrigation Controller has been tested and has demonstrated protection against SAFETY HAZARDS caused by harmful ingress of liquids due to splashing, per IEC 60601-1 and UL 60601-1.

MEDICAL EQUIPMENT
WITH RESPECT TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY
IN ACCORDANCE WITH UL 60601-1 (2003), ANSI/AAMI ES60601-1 (2005),
AND CAN/CSA-C22.2 NO. 60601-1 (2008)
Product contains no serviceable parts.

R_x only

Bard, Davol, SmokEvac, Nezhat-Dorsey, and Endo-Pool are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 2005, 2008, 2011, 2014 C. R. Bard, Inc. All Rights Reserved.

Record the Model Number and Serial Number of the X-STREAM™ Laparoscopic Irrigation Controller and retain for future reference.

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Copyright © 2005, 2008, 2011, 2014 C. R. Bard, Inc. Tous droits réservés.
Filiale de C. R. Bard, Inc.